

Combinatorial Labeling Kits for Long-range nOes Detection

User-friendly solutions for the simultaneous labeling of any combinations of Ala, Ile, Leu, Met, Thr & Val methyl groups.

NMR-Bio kits include the regio-specific and/or stereo-specific labeling of Ile, Leu, and Val residues. Kits are provided with precise protocols extensively tested *in-vivo* to ensure optimal incorporation of isotopes in targeted methyl groups without detectable scrambling in other positions.



2D Methyl-TROSY and 3D ¹³C-edited NOESY spectra of MSG 82 kDa labeled using NMR-Bio kits QLAM-A^βl^{δ1}LV^{proS} & TLAM-I^{δ1}M^βT^γ



Detection of long-range nOes between methyl probes distant by up to 10 Å in large proteins and complexes

NMR-Bio scrambling free labeling solutions are optimized for the extraction of precise and long-range nOe distance restraints between methyl probes in perdeuterated proteins. Compared to standard 2-keto acids, NMR-Bio acetolactate precursors increase sensitivity by up to a factor of 4, allowing for the detection of structurally meaningful long-range and intermolecular nOes restraints.

 $\label{eq:commended kits:} Recommended kits: $$ TLAM-A^{\beta}I^{\delta 1}M^e/I^{\delta 1}LV^{proS}/I^{\delta 1}M^{\epsilon}T^{\gamma}$$ QLAM-A^{\beta}I^{\delta 1}LV^{proS}/A^{\beta}M^{\epsilon}LV^{proS}/I^{\delta 1}LV^{proS}T^{\gamma}/I^{\delta 1}M^{\epsilon}LV^{proS}$$$ PLAM-A^{\beta}I^{\delta 1}M^{\epsilon}LV^{proS}T^{\gamma}$$$ HLAM-A^{\beta}I^{\delta 1}M^{\epsilon}LV^{proS}T^{\gamma}$$$$

References: Kerfah et al., J Biomol NMR. 2015, 61(1):73-82 Kerfah et al., Curr Opin Struct Biol. 2015, 32:113-22 NMR-Bio patented precursors are specifically deuterated and supplied as frozen user-friendly kits ready to use without requirement of any further chemical modification. All NMR-Bio kits are calibrated for addition into deuterated M9 culture medium prior to induction. NMR-Bio kits have been optimized to incorporate ¹³CH₃ or ¹³CHD₂ isotopomers in selected methyl groups of proteins, with the possibility to incorporate a linear ¹³C spin system connecting the specifically labeled methyl



specifically labeled methyl groups to the backbone nuclei. Kits are provided with precise protocols extensively tested *in-vivo* to ensure optimal incorporation of isotopes in targeted methyl groups without detectable scrambling in other positions.

Examples of kits	¹³ CH ₃ groups Labeled	Prices [*] (1 kit for 1 L. of culture)
TLAM-I ^{δ1} LV ^{proS}	Ile ^{δ1} Leu ^{proS} Val ^{proS}	455€
TLAM- $A^{β}I^{\delta 1}M^{ε}$	Ala $^{\beta}$ Ile $^{\delta 1}$ Met $^{\epsilon}$	1385€
ΤLAΜ-Ι ^{δ1} ΜεΤ ^γ	$IIe^{δ1}$ Met ^ε Thr ^γ	1070 €
QLAM-A ^β M ^ε LV ^{proS}	Ala ^{β} Met ^{ϵ} Leu ^{proS} Val ^{proS}	1480 €
QLAM-I ^{δ1} M ^ε LV ^{proS}	$IIe^{\delta 1} Met^{\epsilon} Leu^{proS} Val^{proS}$	865€
$\textbf{QLAM-I}^{\delta 1} \textbf{LV}^{\textbf{proS}} \textbf{T}^{\gamma}$	$IIe^{\delta 1}$ Leu/Val ^{proS} Thr ^{γ}	1005 €
$\textbf{QLAM-A}^{\beta}\textbf{I}^{\delta1}\textbf{L}\textbf{V}^{\textbf{proS}}$	Ala ^{β} Ile ^{δ1} Leu ^{proS} Val ^{proS}	1070 €
$\textbf{PLAM-I}^{\delta 1}\textbf{M}^{\epsilon}\textbf{L}\textbf{V}^{\textbf{proS}}\textbf{T}^{\gamma}$	$IIe^{\delta 1} Met^{\epsilon} Leu/Val^{proS} Thr^{\gamma}$	1415€
$\textbf{PLAM-}A^{\beta}\textbf{I}^{\delta1}\textbf{M}^{\epsilon}\textbf{L}\textbf{V}^{\textbf{proS}}$	Ala $^{\beta}$ Ile $^{\delta 1}$ Met $^{\epsilon}$ Leu/Val ^{proS}	1480 €
$\textbf{PLAM-}A^{\beta}I^{\gamma2}M^{\epsilon}LV^{\text{proS}}$	Ala ^{β} Ile ^{γ2} Met ^{ε} Leu/Val ^{proS}	1505 €
$\textbf{HLAM-A}^{\beta}\textbf{I}^{\delta1}\textbf{M}^{\epsilon}\textbf{L}\textbf{V}^{\textbf{proS}}\textbf{T}^{\gamma}$	Ala ^{β} lle ^{δ1} Met ^{ϵ} L/V ^{proS} Thr ^{γ}	2030€

For any kit including Val-proR & Leu-proR, please inquire !

The listed prices exclude shipping fees and importation tax. Discounts apply for larger quantities. For specific quote contact us at **sales@nmr-bio.com**

www.nmr-bio.com